

California Department of Food and Agriculture

Division of Inspection Services
Center for Analytical Chemistry
3292 Meadowview Rd
Sacramento, CA 95832

Environmental Analysis Section

Jane White- Chemist
Paul Lee- Staff Chemist
Stephen Siegel- Supervising Chemist

Pyrethroids in sediment

- Adapted from: A sonication extraction method for the analysis of pyrethroids, SIU, M.J. Lydy
- Sample may be frozen and batched
- Thaw sample-decant excess water
- Weigh out 20grams of sediment
- Add 5 grams Copper
- Add MgSO_4 mix until consistency of sand
- Add 75mL's 1:1 acetone:hexanes

Pyrethroids in sediment

- Shake sample at 185 rpm for 15 minutes
- Decant through MgSO_4 into boiling flask
- Repeat shaking with 75 additional mL's
- Decant through MgSO_4
- Rinse with 1:1 acetone:hexanes
- Rotovap to approx. 5mL's
- Ready for clean-up procedure

Pyrethroids in sediment water

- Weigh sample plus bottle
- Transfer water portion to separatory funnel
- Add 60mL's hexanes to sample bottle
- Shake for 30 seconds
- Transfer remaining sample to funnel
- Shake funnel for 2 minutes

Pyrethroids in sediment water

- Allow layers to separate
- Drain off lower water layer
- Transfer hexanes through Na_2SO_4 to dry
- Transfer water back into funnel
- Repeat extraction 2 more times
- Collect all hexanes in 500mL boiling flask

Pyrethroids in sediment water

- Rotovap to approx. 5mL's
- Ready for clean-up procedure

Clean-up procedures

- Condition a 2 gram florisil SPE cartridge
- Use 10mL's 15% ether in hexanes
- Follow with 20mL's hexanes
- Load entire extracted sample on SPE cartridge
- Rinse boiling flasks and load onto cartridge
- Elute pyrethroids using 30mL's 15% ether in hexanes

Clean-up procedures

- Evaporate sample to just dryness under gentle nitrogen stream
- Bring to final volume in hexanes
- GC/MS analysis final volume 1.0mL
- GC/ECD analysis final volume 2.0mL

LOD's and LOQ's in sediments

	LOD	LOQ
• Bifenthrin	0.11ppb	1.0ppb
• Fenopropathrin	0.11ppb	1.0ppb
• Permethrins	0.14ppb	1.0ppb
• Fenvalerates	0.11ppb	1.0ppb
• λ-cyhalothrin	0.12ppb	1.0ppb
• λ-cyhalothrin epimer	0.12ppb	1.0ppb

LOD's and LOQ's in sediments

	LOD	LOQ
• Cypermethrins	0.11ppb	1.0ppb
• Deltamethrin	0.10ppb	1.0ppb
• Resmethrin	0.90ppb	1.5ppb

LOD's and LOQ's in water

	LOD	LOQ
• Bifenthrin	1.8ppt	5.0ppt
• Fenopropathrin	1.5ppt	15ppt
• Permethrins	7.7ppt	15ppt
• Fenvalerates	1.8ppt	15ppt
• λ-cyhalothrin	1.2ppt	15ppt
• λ-cyhalothrin epimer	1.1ppt	15ppt

LOD's and LOQ's in water

	LOD	LOQ
• Cypermethrins	1.8ppt	15ppt
• Deltamethrin	1.9ppt	15ppt
• Resmethrin	3.8ppt	15ppt

Stability of pyrethroids in sediment water

- Adding 10mL's hexanes to 1.0 liter sediment water sample prolonged the sample storage stability

Stability of permethrins in sediment water (100 ppt spike) Recovery after 28 days

<u>Compound</u>	<u>w/o hexane</u>	<u>w/hexane</u>
Bifenthrin	54.8%	83.4%
Fenprothrin	77.6%	95.5%
Permethrin-t	14.6%	93.8%
Permethrin-c	42.9%	91.2%
Fenvalerates	44.5%	81.5%

Stability of permethrins in sediment water (100 ppt spike) Recovery after 28 days

<u>Compound</u>	<u>w/o hexane</u>	<u>w/hexane</u>
λ-cyhalothrin	52.8%	86.8%
λ-cyhalothrin		
Epimer	54.2%	90.0%
Cyfluthrin	42.7%	86.6%
Cypermethrin	44.8%	82.3%

Positive attributes of the sediment method

- Simple sample preparation procedure
- Sample preparation uses low cost equipment and commercially available SPE cartridges
- Analysis is performed using GC configured with an ECD detector
- Confirmation is performed using GC configured with a mass selective detector (MSD)

Negative attributes of the sediment method

- Extraction may not be as complete as if done by accelerated solvent extraction (ASE)
- Cleanup procedure may not remove enough matrix contaminants for easy identification on GC/ECD
- May be necessary to raise the reporting limit on "dirty" sediment samples
- Detection limit for some compounds may be too high on the GC/MSD even in SIM mode

Pyrethroid methods

Questions?

CDFG would like to thank:

CALFED

USGS

For allowing us to be part of the pyrethroid method
development