

## USGS Methods: Water, Sediment and Colloids

Michelle Hladik  
USGS California Water Science Center  
Sacramento, CA

## Overview of Methods

- Water
  - Separation of fractions
  - Dissolved
  - Suspended sediments
  - Sorption to containers
- Sediments
  - Bed and suspended sediments
- Colloids
  - SPME




## Pyrethroids Included

- Allethrin
- Bifenthrin
- Cyfluthrin
- Cyhalothrin
- Cypermethrin
- Deltamethrin
- Esfenvalerate
- Fenpropathrin
- Fluvalinate
- Permethrin
- Resmethrin
- Sumithrin (Phenothrin)
- Tetramethrin




## Separation of Sample

Sample large volume of water (100-1000 liters)



## Water Method

- 1-L sample
- Filter (0.7  $\mu$ m glass fiber)
- Solid Phase Extraction (HLB cartridge)
- Bottle rinse (sorbed pyrethroids)
- Extract filter paper (suspended sediments)
- GC/MS and GC/MS/MS (ion-trap)
- Recovery 91-98%
- MDLs 2-5 ng/L




## Pyrethroid Sorption

- Pyrethroids in water sorb to sampling containers (glass or plastic)
  - Up to 50%
  - Analytical and toxicological importance
- Composition of sample influences extent of sorption
  - DOC
  - Suspended sediments
  - Relative surface area of container




## Container Sorption- Varying Removal



- Filtered water, 6 mg/L DOC,
- Equilibrated for 24 hours
- Spiked at 200 ng/L



## Water Mass Balance




- Unfiltered water, 5 mg/L DOC, 13.9 mg/L SSC,
- Equilibrated for 3 days
- Spiked at 400 ng/L



## Pyrethroid Sorption to Containers

- Not all pyrethroids sorb to the same degree
- Allethrin and tetramethrin sorb less  
→ lack similar functional group
- Sorption seems to be "loose association"



## Sorption- Additional Work

- Sorption addressed for analytical purposes not for toxicity or sampling
  - Continue exploring ways to mitigate or address pyrethroid sorption
  - Need concentration organisms are exposed to
- EPA funding to develop SOP for water sampling (FY07-08)
  - Cross-sections, composite samples
  - Autosamplers



## Pyrethroid Water Method- SPE Storage

- SPE cartridge storage
  - Dried
  - Put in freezer
- Tested HLB and C8
  - Found not degradation on either cartridge after 1 month
- Can store cartridges



## Sediment Method

- 5-10 g samples
- MASE extraction at 120 °C, 50 % moisture with DCM:Methanol (9:1)
- Matrix clean-up: Carbon/Alumina stacked cartridges
  - Eluted with DCM
- Sulfur clean-up: GPC or activated Cu
- GC/MS or MS/MS
- Recovery 80-93%
- MDLs 1-5 ng/g



## Sediment Method Development

- **MASE temperature**
  - Did not achieve complete extraction with 100 °C
- **MASE solvent**
  - Started with 1:1 DCM:Acetone
  - More efficient (less matrix) with DCM:Methanol
- **Clean-up**
  - Tried 10 g Florisil
  - No significant reduction of matrix and more time consuming
  - Also took other current use pesticides into consideration




## GC/MS

- **Instrument**
  - Varian
    - 3800 GC
    - Saturn 2000 MS (ion-trap)
- **Column**
  - DB-5MS
    - 30 m × 0.25 mm × 0.25 µm
- **Individual isomer peaks added; reported as total compound**
- **GC/MS/MS**
  - Increases confirmation by decreasing background
  - Does not significantly reduce MDLs



## GC-MS/MS




## MDL Limitations

- **Instrument (ion-trap) is at sensitivity limit**
- **Sediment samples have more matrix problems**
  - Greater concern for multiple isomers
- **Cannot increase sample size (suspended sediments)**




## Colloids- Separation

Sample large volume of water (100-1000 liters)



## SPME

- **Solid-phase microextraction**
- **Fused silica fiber coated with appropriate stationary phase**
- **Used to measure total and bioavailable compounds**



## Colloids- Measurement of Pyrethroids

- Several types of SPME
  - Headspace, liquid
- Negligible depletion SPME
  - Freely dissolved (or bioavailable) fractions
  - Only if equilibrium is not disturbed (EE of less than 5-10%)
  - Sample matrix should not disturb sorption kinetics



## SPME- Parameters

- 20 mL of water
- Fiber
  - Tested 7 and 100  $\mu\text{m}$  PDMS phase
  - 1 cm
- Total pyrethroids
  - Sorption onto fiber at 90 °C for 30 min with agitation
  - Desorption at 275 °C for 3 min
- Negligible depletion/ Bioavailable
  - Based on Liu et al., 2004, ET&C, 23, p 7-11.
  - Sorption at ambient temperature for 10 min with agitation
  - Desorption at 275 °C for 3 min

